

EXALT™ Model B Single-Use Bronchoscope

Suction and best-in-class
visualization with the feel
of a reusable.





Single-use bronchoscopy redefined:

Designed and optimized for a wide array of bedside procedures, EXALT™ Model B features suction capabilities, high-quality imaging, and the feel of a reusable, all with the benefits of single-use.

Evolving single-use technology to elevate patient care:

The EXALT Model B Single-use Bronchoscope is yet another example of our long-held commitment to advancing patient care by evaluating, reimaging, and enhancing single-use imaging technology.



Suction Capabilities

EXALT™ Model B was specifically designed to provide physicians with suction capabilities for use in clearing airways.



High Quality Imaging

Bright LEDs and dynamic light control designed for image optimization.



Familiar Functionality

Placement of controls and working channel are designed for a familiar feel to that of a reusable scope



Sterile, Single Use

EXALT Model B is sterilized and intended for single-use, eliminating the risk of cross-contamination between patients due to ineffective reprocessing.



Robust Suction Performance

Designed with patients, procedures, and suction in mind, EXALT™ Model B is available in three sizes- slim, regular, and large.

- **Slim scope** - the small diameter with suction allows for use within a double-lumen tube during lung isolation procedures
- **Regular scope** - suction capabilities in a 5.0mm scope size, supporting most bedside procedures
- **Large scope** - designed for treatment of patients with voluminous or difficult secretions or that require aspiration of blood clots

[About Model B](#)

[Product Details](#)

[Ordering Information](#)





Familiar Design

Placement of controls and working channel are designed for a familiar feel of a reusable scope

Navigation to target anatomy

The navigation and steerability of EXALT™ Model B was preferred over other single-use bronchoscopes*

Best in Class Visualization

EXALT Model B features dynamic light control designed for optimized, high-quality visualization.

Integrated image capture with Endoscopy Report Writer integration

Capture images during a procedure right from the scope

31,3cm high resolution (2736 x 1824) portable monitor with a 16,5cm live image size, designed for ease of viewing during procedures

*In product evaluations, 15 physicians evaluated EXALT Model B and rated navigation as “clinically preferred” over other commercially available single-use bronchoscopes.

Note: image shown on monitor was captured in a lab setting and not representative of in-human use.



Product Specifications

	EXALT Model B Slim	EXALT Model B Regular	EXALT Model B Large
Field of view	90 °	90 °	90 °
Articulation Angle	180 ° up, 180° down	180 ° up, 180° down	180 ° up, 180° down
Insertion Tube OD	3,8mm	5,0mm	5,8mm
Working Length	60cm	60cm	60cm
Average Working Channel Inner Diameter	1,2mm	2,2mm	2,8mm
Minimum ETT Size Compatibility	5,0mm	6,0mm	7,0mm
Minimum DLT Size Compatibility	35Fr	41Fr	N/A

Ordering Information

Product Code	Description	Unit
M00542760	EXALT Monitor (mount included)	Each
M00542770	EXALT Monitor (mount included) - Alternative languages	Each
M00542790	EXALT Monitor (no mount included) - Alternative languages	Each
M00542780	EXALT Monitor (no mount included)	Each
M00542820	EXALT Monitor mount	Each
M00542751	EXALT Model B single-use bronchoscope – Large	Box 5
M00542731	EXALT Model B single-use bronchoscope – Regular	Box 5
M00542711	EXALT Model B single-use bronchoscope – Slim	Box 5



Power Cords

Product Code	Description	Unit
8500-098BL9	Power Cord (EU) – Plug Type F	Each
8520-098BL2	Power Cord (CH) – Plug Type J	Each
8540-098BL2	Power Cord (DK) – Plug Type K	Each
8550-098BL3	Power Cord (IT) – Plug Type L	Each
8560-098BL3	Power Cord (IL) – Plug Type H	Each
8580-098BL5	Power Cord (ZA) – Plug Type M	Each
9650-098BL7	Power Cord (UK) – Plug Type G	Each



All trademarks are the property of their respective owners. Photos taken by Boston Scientific.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

Caution: The law restricts these devices to sale by or on the order of a physician.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material is not intended for use in France.

Product available in the European Economic Area (EEA) only. Please check availability with your local sales representative or customer service.

**Boston
Scientific**
Advancing science for life™

www.BostonScientific.eu

©2021 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-1008605-AA